

THE BIROBOTICS
INSTITUTE



NeuRo
TouchLab

Sant'Anna
School of Advanced Studies – Pisa

EUCLISA
European University

Hands-on workshop with sensorized mechatronic platforms

Prof. Calogero Oddo, Dr. Fabrizia Auletta

3rd February 2025

The Multisensory Scanner integrates:



Tactile data for
biomechanics detection




Vision data for surface 3D
reconstruction



Acoustic data for bio-material
ultrasonic characterization



AI technologies for data integration
and processing



Our core technology integrates:



Tactile data from uniaxial load cell




Vision data from n.2 RGB-D cameras




Acoustic data from 16 MHz needle ultrasonic probe



AI technologies run on a dedicated workstation (e.g., PC, Desktop)




The Multisensory Scanner ideal workflow




The Multisensory Scanner hands-on workflow



- Object detection in camera reference frame and transformation in scanner reference frame
- Definition of the point to inspect






- Communication with platform
 - Get platform state and data
 - Send point to inspect to platform
 - Set platform state:
 - to move actuators to point
 - to touch the object
 - to stop contact and return home




The Multisensory Scanner hands-on workflow



- Object detection in camera reference frame and transformation in scanner reference frame
- Definition of the point to inspect




The Multisensory Scanner hands-on workflow



Multisensory Scanner research results

Experimental results of the first prototype of the platform proved its ability to localize areas of different biomechanical characteristics inserted in

laboratory-made Agar phantoms, based on mechanical (a, top panel) and acoustic (a, bottom panel) measurements [1], excised animal liver, based on the combined evaluation of amplitude (CIA) and shape (CIS) of the acoustic signal measured by the ultrasound probe [2].






[1] Massari et al. "A mechatronic platform for computer aided detection of nodules in anatomopathological analyses via stiffness and ultrasound measurements." Sensors 19.11 (2019)



[2] Bulletti et al. "An improved strategy for detection and localization of nodules in liver tissues by a 16 MHz needle ultrasonic probe mounted on a robotic platform." Sensors 20.4 (2020)

Multisensory Scanner research results

- The first apparatus was coupled with a novel tele-palpation apparatus [1] to enable the detection of nodules with various distinct stiffness buried in an ad-hoc polymeric phantom [2]
- The effectiveness was proved under two experimental conditions of real-time telepresence:
 - with the platform placed in the visible range of a user
 - with the platform and the user being 50 km apart.



Available data



- Camera images uploaded in shared folder ([link](#))
- Remaining data, shared over UDP connection
 - port 8190
 - Semicolon (;) separated string, containing either 16 or 4017 numbers (double)

Available data – interpretation

Scanner State

Scanner position

Load cell data

Ultrasound data

Timestamp

- Said **S(i)** the received string via UDP:
 - S(0) == 1 if platform is operative, 0 otherwise
 - S(1) == 1 if platform is in the set position, 0 otherwise
 - S(2) == 1 if linear actuator is touching the sample, 0 otherwise
 - S(3) == 1 if linear actuator exceeded its travel range, 0 otherwise
 - S(8) == 1 if linear actuator is in position (0,0,0), 0 otherwise
 - S(9) state id of the PLC state machine
 - S(14) timestamp in seconds of scanner data acquisition
 - S(15) timestamp in seconds of ultrasound data acquisition

Available data – interpretation

Scanner State

Scanner position



Load cell data

Ultrasound data

Timestamp

- Said **S(i)** the received string via UDP:
 - $S(4) + S(11)/1000$ millimeters, is position X
 - $S(5) + S(12)/1000$ millimeters, is position Y
 - $S(6) + S(13)/1000 + S(7)/200$ millimeters, is position Z
 - $S(10) / 27648$ (if $|S(10)| < 27648$) Newton, is Force along Z
 $S(10) / 32767$ (otherwise)
 - $S(16:4017)$ Ultrasound data

To be received data



Scanner next state

Scanner next position



- UDP connection
 - port 8191
- Semicolon (;) separated string


To be received data

Sequence of command to send	Expected behavior
ID *	<ul style="list-style-type: none">• Scanner recognize user and allow user control
ID; S0; X; Y; Z **	<ul style="list-style-type: none">• Scanner moves to position X, Y, Z
ID; S1	<ul style="list-style-type: none">• Scanner start moving downward and stop if in contact
ID; S2	<ul style="list-style-type: none">• Scanner stop contact and returns home

* Each user will set their own ID

** X, Y, Z represent the position in millimeters

Multi sensory scanner – state machine



Other resources available

- Transformation matrix in
 - folder name “test_cal”
- Images and pointcloud of the empty platform in
 - folder name “test_0_empty_20250203” for rear RGBD camera
 - folder name “test_1_empty_20250203” for front RGBD camera
- Remember:
 - D405** identifies rear camera
 - D435** identifies front camera

LET'S TRY IT

A word cloud in the shape of a spiral, containing various terms related to haptics, touch, and prosthetics. The words are arranged in a clockwise spiral pattern, with larger words in the center and smaller words towards the outer edges. The terms include:

- Artificial
- Fingertip
- Skin
- Medical
- robotics
- Haptics
- Neuroprosthetics
- Restoration
- MEMS
- Wearable sensors
- Industry 4.0
- Telepresence
- Haptic displays
- Heart biomechanical sensing
- Microneurostimulation
- Microneurography
- Human-machine cwork
- Sensory substitution
- Tactile actuators
- collaborative robotics
- Sensory augmentation
- FBGs
- consumer electronics
- Mechanoreceptors
- Wearable technologies
- Cardiorespiratory monitoring

Object detection suggestions

- Read information about the scene, including camera settings, images, and point cloud data.
- Load calibration files to understand how cameras are positioned.
- Load a base reference mask that represents the platform.
- Verify that the cameras used for capturing the scene are the same as those in the calibration.
- Ensure they have the same type and serial number to avoid mismatches.
- For each captured scene, get the point cloud data (3D points of the object).
- Use stored transformation data to align it to the platform.
- Apply the necessary transformations so that object points are placed in the platform's reference frame.
- Compare each captured image with a reference image from calibration.
- Identify the region where the object is located.
- Extract a mask (a filtered version of the point cloud) that contains only the object.

Available data - recap

Scanner State

Scanner position

Load cell data

Ultrasound data

Timestamp

tri-axial motorized platform	isReady_Festo	1 if platform is operative, 0 otherwise	bit
	isArrived_Festo	1 if platform is in the set position, 0 otherwise	bit
	isTouching_SMAC	1 if linear actuator is touching the sample, 0 otherwise	bit
	isMaxStroke_SMAC	1 if linear actuator exceeded its travel range, 0 otherwise	bit
	posx	x-axis position (integer part)	mm
	posy	y-axis position (integer part)	mm
	posz	z-axis position (integer part)	mm
	posz2	voice coil z-axis position	a.u.
	isHome_SMAC	1 if platform is in position (0,0,0), 0 otherwise	bit
	state_Festo	state id of the PLC state machine	number
	Fz	normal force	a.u.
	posx	x-axis position (fractional part)	mm
	posy	y-axis position (fractional part)	mm
	posz	z-axis position (fractional part)	mm
Ultrasound	CPXEts	timestamp of data collection	seconds
	USts	timestamp of data collection	seconds
	0-3999	ultrasound signal	a.u.

